AP BIOLOGY DIVERSITY OF LIFE ACTIVITY #2

NAME			

DATE_____HOUR____

PROKARYOTES

PART I: GRAM STAINING OF BACTERIA

INTRODUCTION

Gram stain is commonly used to assist in bacterial identification. This stain, first developed in 1884, separates bacteria into groups, depending on their reaction to this stain. Bacteria react by testing **gram-positive**, **gram-negative**, **gram-variable**, with the first two groups being the most common. The response of cells to the stain is due to differences in their cell walls. Cell walls of gram-positive bacteria are predominantly a complex protein-carbohydrate polymer, **peptidoglycan**, located outside the plasma membrane. Cell walls of gram-negative bacteria are thinner, with less peptidoglycan present. The wall in gram-negative bacteria is a polysaccharide-protein-lipid complex. Studies of bacterial taxonomy have show that these differences define major taxonomic groups.

Gram stain relies on the use of two dyes, crystal violet (purple) and safranin (pink). Cells with the peptidoglycan wall retain the purple dye and are **gram-positive**. Those cells that retain the pink safranin are **gram-negative**.

In this lab, you will prepare and stain slides of five different bacterial species to determine how each reacts to gram stain.

MATERIALS

<u>At the "smear" station:</u>

Bunsen Burner Burner striker Inoculating loop Beaker of distilled water Bacteria cultures Clean microscope slide

At staining station:

Staining jars of: Crystal Violet Gram's Iodine 95% Ethyl Alcohol Safranin 1000 mL beaker with tap water

PROCEDURE -- "SMEAR" PREPARATION

In this section of the lab, your group will prepare a smear of each of the bacteria cultures. Make sure you carefully follow the steps below to avoid contamination of the bacteria cultures and of the lab area.

1. Carefully light the burner. Adjust the burner so you have a blue flame.

2. Sterilize the end of the inoculating loop by placing it in the flame of the burner until the wire turns bright red.

CAUTION: Use care with the open flame of the burner. Do not touch the hot end of the inoculating loop. The loop may still be quite hot after the red glow has faded.

- 3. With a wax pencil, label one end of the slide #1. Place a loopful of distilled water on a clean, dry microscope slide. Sterilize the loop. Allow the loop to air cool before proceeding to the next step.
- 4. Lift the petri dish lid of the bacteria culture. Gently touch the sterile inoculating loop to a bacteria colony. The bacteria should be barely visible on the loop. Replace the petri dish lid. Spread the bacteria sample on the microscope slide by mixing it with the drop of water.
- 5. Sterilize the loop. Allow the loop to air cool before proceeding to the next step.
- 6. Allow the smear to dry by waving the slide over the Bunsen burner flame. To prevent the smear from cracking and peeling from the slide, do not heat the slide more than the fingers can comfortably stand.
- 7. After the film is thoroughly dry, FIX it by passing it quickly through the flame two or three times. Be sure the specimen side of the slide is up when fixing. Allow the slide to air cool before proceeding to the staining steps. Turn the burner off.
- 8. Repeat steps 1 through 7 for each culture of bacteria. Be sure to label each slide with the number of the bacteria culture.

PROCEDURE -- GRAM STAIN


- 9. Use a clothes pin to hold the slide during the staining procedure. Place the smear in the crystal violet for 60 seconds.
- 10. Rinse the slide quickly (a few seconds) in tap water and immediately place it in the Gram's iodine solution for 60 seconds.
- 11. Rinse the smear gently in tap water and dry it by blotting with a clean paper towel. (**CAUTION:** Blot gently so that the bacteria are not rubbed from the slide.)
- 12. Place the smear in 95% ethyl alcohol, a decolorizing agent, for 30 seconds. Blot the slide again with a clean paper towel.
- 13. Immerse the slide in safranin, a couterstain, for ten seconds.

- 14. Wash the slide in tap water and blot it dry.
- 15. Examine the slide under high power of the microscope. The gram positive bacteria cells will appear violet in color and the gram negative cells will be red or pink in color.

OBSERVATIONS

16. Complete the following chart. Use the pictures at the bottom of the page to help you identify the cell shape and cell arrangement.

Culture	1	2	3	4
Drawing (in color)				
Magnification				
Cell Size				
Cell Shape				
Cell Arrangement				



17. Use the information below and the information you collected in this lab to match the bacterial culture (1, 2, 3, or 4) to the correct species.


Bacillus cereus (gram-positive, rod)
 Micrococcus luteus (gram-positive, coccus)
 Rhodospirillum rubrum (gram-negative, spiral)
 Serratia marcescens (gram-negative, rod)

PART II: QUESTIONS

1. At present, the placement of all prokaryotes into Kingdom Monera is being challenged. Why?

2. The shape and arrangement of bacterial cells are used to identify bacteria. On page 5 are drawings of several disease-causing bacteria. Use the key below to identify the bacteria in each picture.

1	a. Bacillus bacteriab. Coccus bacteriac. Spirillum bacteria	go to 2 go to 4 go to 7
2.	a. Bacteria arranged in chainsb. Bacteria appear singly	Rat Bite Fever Stain A go to 3
3.	a. Bacteria with swollen endosporeb. Bacteria without swollen endospore	Tetanus Anthrax
4.	a. Bacteria appear in pairsb. Bacteria appear in chainsc. Bacteria appear in clumps	go to 5 Strep Throat go to 6
5.	a. Bacteria have pointed endsb. Bacteria without pointed ends	Pneumonia Gonorrhea
6.	a. Bacteria appear in groups of 4 (Sarcina)b. Bacteria appear in clusters	Human skin bacteria Staph infection
7.	a. Bacteria with flagellab. Bacteria without flagella (Spirochete)	Rat Bite Fever Stain B Syphilis



- 3. What are three functions of the cell wall in prokaryotes?
- 4. What are peptidoglycans?
- 5. How do antibiotics affect peptidoglycans?

6. Identify each of the following as true of gram-positive (+) or gram-negative (-) prokaryotes.

Simpler cell walls
Cell walls contain large amounts of peptidoglycans
Cell walls more complex
Cell walls contain smaller amounts of peptidoglycans
Have a membrane containing lipopolysaccharides outside the cell wall Stain purple with Grams stain
Stain pink with Grams stain
Why are gram-negative bacteria more often disease causing than gram- positive?

8. List the 3 types of motility used by motile prokaryotes.

9. Match the term with the correct definition

7.

Movement toward or away from a stimulus	A. Chemotaxis
Movement toward or away from	B. Phototaxis
Movement toward or away from a certain chemical	C. Taxis

10. Explain why a photosynthetic bacterium would exhibit a positive phototaxis. A positive taxis occurs when the bacterium moves toward the stimulus.

The prokaryotic flagellum is not homologous to the eukaryotic flagellum. Use 11. the key below to indicate if the characteristics listed below are true of prokaryotic or eukaryotic flagella.

P = Prokaryotic	E = Eukaryotic	
smaller width	larger width	
covered with plasma	not covered with	ano
composed of flagellin	•	
9+2 arrangement of	no 9+2 arrangemer	ıt
Rotating motion	Undulating motion	
Filament is solid	Made of microtubule	es
ch the term with the correct de	escription of definition.	
A. Binary fission B. Endospores	C. Genophore D. Plasmids	
Prokaryotic chromosome		
Usually double-stranded;	forms a ring; found in nucleoid region	
Small rings of DNA		
Consist of only a small null	mber of genes	
Endows cell with antibiotic unusual nutrients	resistances or the ability to metabolize	ž
Method of reproduction in	prokaryotes	
Pesistant cells that contain by a durable wall	n one copy of the chromosome surroun	ded
Allows some prokaryotes t	to withstand harsh conditions	
	 smaller width covered with plasma membrane composed of flagellin 9+2 arrangement of microtubules Rotating motion Filament is solid the term with the correct detection of the term with the correct detection of term with the correct detection of term with the term with the term with the correct detection of term with the correct detection of term with the term with the term with the term with the correct detection of term with the term with term with the term with term with term with the term with	smaller width

- 14. You have a bacterial infection and the doctor gives you an antibiotic that blocks the protein synthesis of the bacteria in your body. What effect does this antibiotic have on the protein synthesis in your cells? Explain.
- 15. Match the definition with the correct term. _ Change in DNA base sequence A. Conjugation Origin of a new trait B. Mutation _____ Major source of variation in C. Transduction prokaryotes _____ Direct transfer of genes between D. Transformation two prokaryotes _____ Cell takes in and incorporates DNA from the environment Transfer of genes between cells via a virus 16. Use the key below to identify the nutrition type being described. PA = Photoautotrophs **PH = Photoheterotrophs** CA = Chemoautotrophs CH = Chemoheterotrophs _____ Energy source is light; Carbon source is organic Energy source is oxidation of inorganic substances (H₂S, NH₃, Fe²⁺); Carbon source is CO₂ Energy source is organic; Carbon source is organic Energy source is light; Carbon source is CO₂ Two modes of nutrition restricted to certain prokaryotes; not found in other organisms Many prokaryotes, many protists, all fungi, and all animals

Most plants

- 17. Identify each of the following as true of **S**aprobes or **P**arasites
 - _____ Decomposers
 - _____ Absorb nutrients from dead organic matter
 - _____ Absorb nutrients from the body fluids of a living host
 - _____ The Chlamydias live inside the cells of a host; leading cause of blindness and causes the most common form of STD in the US.
 ____ Prokaryotes found in or on food
 - _____ Chemoheterotrophs
- 18. Match the description or definition with the correct term.

Poisoned by oxygen	A. Facultative anaerobes
Cannot grow without oxygen	B. Obligate aerobes
Cannot use oxygen	C. Obligate anaerobes
Use oxygen if it is present; uses	

- fermentation oxygen is absent
- 19. Listed below are characteristics of each group of Domain Archaea. Use the key below to indicate which group the characteristic applies.

Μ	=	Methanogens
EH	=	Extreme halophiles
Т	=	Thermoacidophiles

- _____ Found in mineral springs and thermal volcanic vents on ocean floor
- _____ Found in waters of extreme salinity (15 to 20%)
- _____ Some are important symbionts in termites & other herbivores
- _____ Some are important decomposers in marshes and swamps
- _____ Contain the pigment bacteriorhodopsin
- _____ Use H_2 to reduce CO_2 to CH_4
- _____ Habitat is hot (60°C to 80°C) and acidic (pH 2 to 4)

How i	s photosynthesis in cyanobacteria similar	to that in plants?
	ria are important in the chemical cycling t acterial group with the correct role.	hat occurs on earth. Matc
	Convert N ₂ to NH ₃ and NO ₃	A. Autotrophic bacter
	$_$ Convert NH ₃ to NO ₂ ⁻ & NO ₃	B. Decomposers
	Convert NO_3^- to N_2	C. Denitrifying bacter
	Fix CO ₂	D. Nitrifying bacteria
	Support food chains	E. Nitrogen-fixing bacteria
	Some produce oxygen gas	Dacteria
	Break down dead organisms and waste	2
	Return elements (C and H) to the envir	ronment
Match	the definition or description with the cor	rect term.
	Ecological relationship between organisms of different species that are in direct contact	A. Commensalism
	Organisms involved in symbiosis	B. Host
	Usually the larger symbiont	C. Mutualism
	Symbiotic relationship where both symbionts benefit from the relationship	D. Parasitism
	Symbiotic relationship where one symbiont benefits and the other is neither harmed nor helped	E. Symbiont
	Symbiotic relationship where one symbiont benefits at the expense of the other	F. Symbiosis

- 23. Legumes are a group of plants that contain nitrogen-fixing bacteria in root nodules.
 - a. What does each organism get from this relationship:

The plant gets:	
• •	

The bacteria get:_____

- b. What type of symbiotic relationship (mutualism, commensalism, or parasitism) does this represent?
- 24. Humans exploit bacteria and other prokaryotes for scientific and commercial purposes. List 5 of those uses.